

Mouse Direct PCR Kit

Description

The Mouse Direct PCR Kit provides a fast preparation and PCR amplification that is specifically designed for mouse genotyping. The Buffer L and Protease Plus rapidly digest mouse **tail**, **ear** and **toe** to release intact genomic DNA that can be used directly as the template for PCR amplification. By using this kit, the digestion process only takes **15 minutes**. In addition, the 2x M-PCR OPTI Mix (which includes an optimized Taq Polymerase) ensures high amplification efficiency of target DNA.

Components

| Contents | Cat#:B40013 | Cat#:B40015 |
|--------------------|-------------|-------------|
| Buffer L | 20 ml | 25 mL x 2 |
| Protease Plus | 0.4 mL | 1 mL |
| 2 x M-PCR OPTI Mix | 1 mL x 2 | 1 mL x 5 |

- Buffer L: Lysis buffer
- Protease Plus: For rapid and efficient digestion of mouse tissue in only 15 minutes!
- 2x M-PCR OPTI Mix: Includes Bimake's optimized Taq DNA polymerase, dNTPs, MgCl₂, and reaction buffer.


Storage

- 2x M-PCR OPTI Mix and Protease Plus should be stored at -20°C for 2 years. If the PCR Mix is to be used frequently, it can be stored at 4°C for up to 10 days.
- Buffer L should be stored at 4°C.

Notice


1. Ensure all tools are clean to **avoid cross-contamination** when collecting samples.
2. All reagents have been optimized for **use together** and any modifications or alternative uses are prohibited.
3. Before each step, check to make sure every reagent is fully **liquid and resuspended** prior to use.
4. When performing a digestion, the cut fragment must be **completely immersed** in the digestion solution.

Protocol




Add **protease mixture** to the tissue tub

..... **Tissue Digestion**



Incubate at **55°C** for **15** minutes to release genomic DNA

..... **Protease Inactivation**



After digestion, incubate at **95°C** for **5** minutes to inactivate the proteases

..... **PCR**



After centrifuge, 12,000rpm for 10 min, the supernatant of digested solution can be used as the template for PCR



Order & Inquiry

Tel: (713)732-2181 Fax: +1-866-747-4781
E-mail: order@bimake.com



Order & Inquiry

Tel: +49-89-46148500 Fax: +49-89-461485022
E-mail: eu.order@bimake.com

Genomic DNA Preparation

1. Place the mouse **tail, ear, or toe** in a 1.5 mL centrifuge tube.
2. **Thoroughly mix** 100 µL of fresh **Buffer L** with 2 µL of **Protease Plus** for a single sample in a separate tube.
3. Add the **protease mixture** to the mouse tissue tubes **with the tissue cut end submerged in it**, then incubate at **55°C for 15 minutes** (incubation times may vary, depending on samples digestion rates).
4. After the digestion process, incubate at 95°C for 5 minutes to inactivate protease. The tissue lysate can now be used as a PCR template.

PCR Genotyping

5. Add ddH₂O primers, template, and 2 x M-PCR OPTI Mix into a PCR tube according to the recommended concentrations. Give the mixture a quick spin in the centrifuge and load into PCR amplifier to begin amplification.

| PCR Reaction Components | 20 µL Reaction Volume (µL) | 50 µL Reaction Volume (µL) |
|-------------------------|----------------------------|----------------------------|
| ddH ₂ O | 8 | 21 |
| Forward Primer (10 µM) | 0.5 | 1 |
| Reverse Primer (10 µM) | 0.5 | 1 |
| Template | 1 | 2 |
| 2 x M-PCR OPTI Mix | 10 | 25 |

Note: The volume of template is Variable.

| Temperature (°C) | Time | Cycles |
|------------------|-------------------|--------|
| 94 | 5 min | 1 |
| 94 | 20 sec | 35 |
| 50-65 | 30 sec | |
| 72 | X min (2 kb /min) | |
| 72 | 5 min | 1 |
| 12 | -- | 1 |

Trouble Shooting

Please review the following for trouble-shooting options when you encounter technical difficulties. Alternatively, feel free to contact Bimake technical support directly.

| Problem | Potential Cause(s) | Suggestion(s) |
|--|--|--|
| No amplification product in test or control samples | Amplification reaction was incorrectly set up | Optimize the proper reaction set up |
| | Improper storage has led to loss of activity of PCR reagents | Replace all components with fresh reagents |
| | Primers are not optimal and did not anneal | Redesign primers |
| Amplification worked in the control samples, but not in test samples | Digestion was incomplete | Extend digestion time up to 30 minutes at 55°C |
| | Lysis solution was mixed with PCR mixture for too long | Collect fresh mouse tail samples for genomic DNA extraction |
| | The quantity of the amplification product was not sufficient | Increase the number of PCR cycles to 35-40 to yield more amplification product |
| Non-specific amplification product(s) | Annealing temperature was too low | Increase the annealing temperature |
| | The number of PCR cycles was too high | Decrease the number of cycles to 30-35 |
| | Primer concentration was too high | Decrease primer concentration |
| | Template concentration was too high | Dilute template in purified H ₂ O or TE buffer |



Order & Inquiry

Tel: (713)732-2181 Fax: +1-866-747-4781
E-mail: order@bimake.com



Order & Inquiry

Tel: +49-89-46148500 Fax: +49-89-461485022
E-mail: eu.order@bimake.com