


Poly FLAG Peptide lyophilized powder

Column: 4.6mm*250mm,Inertsil ODS-SP
 Solvent A: 0.1%Trifluoroacetic in 100% Acetonirile
 Solvent B: 0.1%Trifluoroacetic in 100% Water
 Gradient:

	A	B
0.01min	12%	88%
25.00min	37%	63%
25.01min	100%	0%
30.00min	Stop	

 Flow rate: 1.0ml/min
 Wavelength: 220nm
 Volume: 10 μ l



Rank	Time	Quantity	Area	Height
1	10.608	0.4272	28499	3514
2	10.792	95.33	6359839	813987
3	10.983	4.244	283171	44728
Total		100	6671509	862229